



Be Right™



Manganese test kit, model MN-5

Produktnr:

146700

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges

!

!

!

!

!

Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. An economical solution, this kit uses the simple, accurate color disc method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	0 - 3.0 mg/L Mn
Metodnamn:	Color disc/Cold periodate
Modell:	MN-5
Parameter:	Manganese
Plattform:	Color Disc
Smallest increment:	0.1
Standardmätområde:	0 - 3.0 mg/L Mn
Test-kit modell:	MN-5