



Acetatbuffertlösning, pH 4,0, 100 mL

Produktnr:

1490932

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Acetate Buffer Solution, pH 4.0. Total and/or Free Chlorine. Method: Amperometric. pk/100 mL.

Specifikationer

Antal tester :	100
Instrument:	AutoCat [™] 9000
Mätområde:	0.0012 - 5.000 mg/L Cl ₂
Metodnamn:	Amperometric Titration
Metodnummer:	8209
Plattform:	Lösning