



HydraVer 2 hydrazinreagens, 500 mL

Produktnr:

179049

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



For hydrazine determination in boiler waters by the p-Dimethylaminobenzaldehyde method. Hach method 8141. Range: 4-600 µg/L N₂H₄. (Also used in visual test kits. Range: 0.02-1.00 mg/L) 500 mL bottle

Specifikationer

Antal tester :	100
Instrument:	All except DR 900, PC II
Lagringsförhållanden:	10 - 25 °C
Mätområde:	4 - 600 µg/L N ₂ H ₄
Metodnamn:	p-Dimethylaminobenzaldehyde
Metodnummer:	8141
Parameter:	Hydrazine
Plattform:	Lösning
Volym/förpackningsstorlek:	500 mL