



Chromium color disc test kit, model CH-8

Produktnr:

183400

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Widely used as a corrosion inhibitor in cooling towers and engine cooling systems, chromium is also used in surface finishing, leather tanning, and other industries. The most widely used form, hexavalent chromium, is very toxic and strict limits are placed on discharges. Utilizing the color disc method to determine chromium (Range: 0.1-1.5 mg/L), this kit contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	0 - 1.3 Cr ⁶⁺
Metodnamn:	Color disc/ Diphenylcarbazide
Modell:	CH-8
Parameter:	Chromium Hexavalent
Plattform:	Color Disc
Smallest increment:	0.1
Test-kit modell:	CH-8