



Be Right™



Hydrazine color disc test kit, model HY-2

Produktnr:

184900

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method with a range of 0.02-1 mg/L, this kit contains 300 tests.

Measures Hydrazine, a common oxygen scavenger in boiler water, using the simple color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	300
Kapslingstyp:	D
Mätområde:	0 - 1 mg/L N ₂ H ₄
Metodnamn:	Color disc/ p-Dimethylamino- benzaldehyde
Modell:	HY-2
Parameter:	Detergents
Plattform:	Color Disc
Smallest increment:	0
Standardmätområde:	0 - 1 mg/L N ₂ H ₄
Test-kit modell:	HY-2