



Bottle, BOD, 60 mL capacity, marked 30 mL graduation line

Produktnr:	190902
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 mL sizes).

Specifikationer

Body Material:	Borosilicate Glass
Kapacitet :	60 mL
Material:	Borosilicate Glass
Package size UOM:	1pk
Volym/förpackningsstorlek:	1 each