



Fenylarsinioxid (PAO) standardlösning, 0,00564 N, 1 L

Produktnr:

199953

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



For chlorine determination by amperometric titration methods. Note: also available in Digital Titrator cartridges. 0.00564 N, 1 liter bottle

Specifikationer

Instrument:	AutoCat [™] 9000
Mätområde:	0.0051 - 5.000 mg/L Cl ₂
Metodnamn:	Amperometric Titration
Metodnummer:	10025
Plattform:	Lösning