



Test kit, acidity, model AC-DT

Produktnr:

2064000

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Acidity is usually caused by the presence of mineral acids, salts of strong acids, or free carbon dioxide. In surface or ground waters, acidity may form natural substances or industrial pollution. Excessive levels cause corrosion and can be detrimental to fish. Utilizing the Hach Digital Titrator, the Model AC-DT kit is used to determine acidity concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-4,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specifikationer

Antal tester :	100
Fotnot :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Kapslingstyp:	F
Mätområde:	10 - 4,000 mg/L CaCO ₃
Metodnamn:	Digital Titrator/Sodium hydroxide
Minsta ökningssteg:	0.1-1.10
Parameter:	Acidity - (P or Total) as CaCO ₃
Plattform:	Digital Titrator
Test-kit modell:	AC-DT