



Cone, Imhoff settling, 1 L, polycarbonate, w/ screw closure

Produktnr:	206700
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Fully graduated cones for determination of small amounts of settleable matter in sewage by Imhoff Sedimentation Method. Glass cone features reinforced tip for ruggedness, while the Polycarbonate cone features a leak-proof, removable closure for withdrawal of sediment for analysis. Support ring, support stand, and cone holders available and sold separately.

Specifikationer

Kapacitet :	1 L
Volym/förpackningsstorlek:	1 each