



Cylinder, Graduated, Mixing, Glass, 10 mL ± 0.1 mL, 0.2 mL divisions, glass stopper #13

Produktnr:	2088638
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

This glass graduated cylinder is calibrated “to contain” with a molded in base. Each includes glass penny head stopper. Graduations are in durable white enamel bake.

Specifikationer

Certifieringar :	Class B
Kapacitet :	10 mL
Material:	Glass
Mätmetod:	To Contain
Typ:	Single Metric, Mixing
Volym/förpackningsstorlek:	1 each