



Cylinder, graduated, mixing, Glass, 50 mL ± 0.4 mL, 1.0 mL divisions, glass stopper #16

Produktnr:	2088641
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

This tall form glass graduated cylinder is calibrated “to contain” with a molded in base. Each includes glass penny head stopper. Graduations are in durable white enamel bake.

Specifikationer

Certifieringar :	Class B
Kapacitet :	50 mL
Material:	Glass
Mätmetod:	To Contain
Typ:	Single Metric, Mixing
Volym/förpackningsstorlek:	1 each