



**Be Right™**



## **Flask, Erlenmeyer, glass w/screw cap, 50 mL**

<b>Produktnr:</b>	<b>2089641</b>
SEK Pris:	Kontakta oss
Ring för leveransdatum	

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

---

### **Specifikationer**

Enhet:	1 each
Kapacitet :	50 mL