



Flask, Erlenmeyer, Polycarbonate Capacity 125 mL

Produktnr:	2089743
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Graduated (Molded-in graduations) flask features: clarity, toughness and strength. Autoclavable; however, repeated autoclaving may result in loss of mechanical strength. #5 Neoprene stoppers (P/N 14808-05; pk/12) available separately. 125 mL.

Specifikationer

Enhet:	1 each
Kapacitet :	125 mL
Material:	Polycarbonate