



Flask, Erlenmeyer, polymethylpentene capacity 250 mL, 4/pk

Produktnr:	2089876
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

Specifikationer

Enhet:	4/pk
Kapacitet :	250 mL
Material:	Polymethylpentene
Typ:	Erlenmeyer