



Buffert, fosfattyp, pH 7, för Amperometric Chlorine, 100 mL MDB

Produktnr: 2155332
SEK Pris: Kontakta oss
Inget leveransdatum anges

For determination of chlorine by amperometric titration. APHA formulation 100 mL bottle

Specifikationer

Antal tester :	100
Instrument:	AutoCat [™] 9000
Mätområde:	0.0051 - 5.000 mg/L Cl ₂
Metodnamn:	Amperometric Titration
Metodnummer:	8334
Plattform:	Lösning