



Test kit, nitrite, nitrogen, model NI-15

Produktnr:

2182000

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Simple, low-cost color disc method is the most accurate visual test kit method.

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	Single Range
Metodnamn:	Color disc/Diazotization
Modell:	NI-11
Parameter:	Nitrite
Plattform:	Color Disc
Smallest increment:	0.01
Standardmätområde:	0 - 0.40 mg/L NO ₂ ⁻ -N
Test-kit modell:	NI-15