



Be Right™



Test kit, formaldehyde, model FM-1

Produktnr:

2183100

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges

Formaldehyde solutions are used for maintenance of reverse osmosis systems and as disinfectants. The Model FM-1 measures higher concentrations by drop count titration. Contains 100 tests.

Utilises drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

| | |
|---------------------|-------------------------------------|
| Antal tester : | 100 |
| Kapslingstyp: | D |
| Mätområde: | 0,05 - 1% & 0,5 - 10% |
| Metodnamn: | Drop count titration/ |
| | Thymolphthalein |
| Parameter: | Formaldehyde - As CH ₂ O |
| Plattform: | Drop Count Titration |
| Smallest increment: | 0.05% |
| | 0.5% |
| Test-kit modell: | FM-1 |