



Reagenssats kiseldioxid, lågt intervall, 25 mL

Produktnr:

2256900

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



For low range silica determination by the Heteropoly Blue method. 25 mL sample Approximately 100 tests

Specifikationer

Antal tester :	100
Mätområde:	0.010 - 1.600 mg/L SiO ₂ and 3 - 1,000 µg/L SiO ₂
Metodnamn:	Heteropoly Blue
Specialfunktioner:	LR, 25 mL