



Rör för destillationsapparat för allmänt ändamål, rund botten, 500 mL

Produktnr:

2264900

SEK Pris:

Kontakta oss

Inget leveransdatum anges

Designed for use with the General Purpose Distillation Apparatus. Round bottom, 500 mL capacity. Long neck with branch and thermometer port. Neck extends 90 mm above flask body and branch extends 150 mm above flask body. Both neck and branch have 20-400 screw-thread finish. Thermometer port has 13-425 screw-thread finish.

Specifikationer

Volym/förpackningsstorlek: 1 each