



Reagenssats, koldioxid, digital titrator

Produktnr:

2272700

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



For determination of Carbon Dioxide by Digital Titrator titration. Hach Method 8205. Range: 10 - 1000 mg/L as CO₂ Set includes: Phenolphthalein Pillows and two Sodium Hydroxide Digital Titrator Cartridges (0.3636 N and 3.3636 N). Approximately 40-100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikationer

Antal tester :	40 - 100
Mätområde:	10 - 1000 mg/L CO ₂
Metodnummer:	8205
Plattform:	Digital Titrator