



Be Right™



Test kit, hydrogenperoxide, model HYP-1

Produktnr:

2291700

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Model HYP-1 Hydrogen Peroxide Test Kit (as H_2O_2) utilizes drop count titration and contains approximately 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	0,2 - 2,0 mg/L & 1 - 10 mg/L
Metodnamn:	Drop count titration/ Thiosulfate
Parameter:	Hydrogen Peroxide - As H_2O_2
Plattform:	Drop Count Titration
Smallest increment:	0.2 1
Test-kit modell:	HYP-1