



Digital titratorpatron, järnhaltigt etylendiammoniumsulfat, 0,00564 N

Produktnr:

2292301

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



For chlorine determination by Digital Titrator titration, using DPD Free or Total Chlorine Reagent indicator. Hach Method 8210. Range: Up to 3.00 mg/L Cl₂. 0.00564 N FEAS. Digits x 0.01 = mg/L Cl₂ in a 25 mL sample.

Specifikationer

Mätområde:	0 - 3.00 mg/L Cl ₂
Metodnamn:	DPD-FEAS
Parameter:	Chlorine, Free
Plattform:	Digital Titrator
Titrant:	0.00564 N FEAS