



Test kit, iron color disc, model IR-21

Produktnr:

2299300

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Iron color disc/TPTZ low range test kit contains 100 tests.

Simple, fast, color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	0 - 2.0 mg/L Fe
Metodnamn:	Color disc/TPTZ
Modell:	IR-21
Parameter:	Iron TPTZ
Plattform:	Color Disc
Smallest increment:	0.01
	0.05
Standardmätområde:	0 - 2.0 mg/L Fe; 0 - 1.2 mg/L Fe
Test-kit modell:	IR-21
Typ:	Double Metric