



Presence/absence test, pk/12

Produktnr:	2323212
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Presence/Absence test for coliforms; Disposable; pk/12.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

Specifikationer

Material:	Plastic
Volym/förpackningsstorlek:	12 /pk