



## Filter holder, for vacuum filtration, glass

|                           |              |
|---------------------------|--------------|
| Produktnr:                | 234000       |
| SEK Pris:                 | Kontakta oss |
| Inget leveransdatum anges |              |

Beaker-style funnel (300 mL) with 47 mm membrane support for filtration of aqueous samples. Fritted borosilicate glass membrane support for even filtration across surface of membrane. Funnel attaches evenly to base via aluminum clamp. Filter flask sold separately.

---

### Specifikationer

|                            |                                                                                                                       |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------|
| Diameter:                  | 47 mm                                                                                                                 |
| Vad ingår i leveransen:    | Borosilicate glass funnel and base; anodized aluminum clamp; support screen (sintered glass); rubber stopper (No. 8). |
| Volym/förpackningsstorlek: | 1                                                                                                                     |

---

### Vad finns i förpackningen

Borosilicate glass funnel and base; anodized aluminum clamp; support screen (sintered glass); rubber stopper (No. 8).