



Be Right™



Test kit, oxygen scavenger, color disc, model OS-1, 50 tests

Produktnr:

2349500

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method, this kit contains 50 tests.

Simple and fast color disc method

Measures common oxygen scavengers: DEHA, hydroquinone, ISA, and MEKO

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	50
Kapslingstyp:	R
Mätområde:	0 - 300 µg/L DEHA (and other Oxygen Scavengers)
Metodnamn:	Color disc/Iron reduction
Modell:	OS-1
Parameter:	Oxygen Scavenger - As DEHA, hydroquinone, erythorbic acid, and methylethylketoxime
Plattform:	Color Disc
Test-kit modell:	OS-1