



Filter, acrodisc, pore size 0.2 µm

Produktnr:	2415600
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

For filtration of dilute biological fluids. Sterilized by gamma irradiation to eliminate potential contamination by EtO residuals. Inlet/Outlet connections: 13/25 mm: Female Luer-Lok inlet, standard male slip luer outlet. Pore size: 0.2 µm. Diameter: 25mm.

Specifikationer

Diameter:	NA
Filter Type:	Filter, acrodisc
Material:	Polysulfone
Packing:	NA
Porstorlek:	0.2 µm