



Reagenssats, klorid, titrering med droppräkning

Produktnr:

2437600

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



For determination of chloride by argentometric drop-count titration. Replacement reagent set for 8-P Test Kit (144001). Set includes: Chloride 2 Indicator Powder Pillows and Silver Nitrate Titrant (0.0493 N titrant). Approx. 100 tests.

Specifikationer

Mätområde:	20 - 400 mg/L Cl ⁻
Metodnamn:	Silver Nitrate Titration