



## Reagenssats, klor (totalt), CN-65 testsats

**Produktnr:**

**2438600**

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



For determination of Total Chlorine by Drop Count Titration. Replacement reagent set for CN-65 Test Kit (225401). Set includes: Sulfamic Acid Powder Pillows, Sulfite 1 Reagent Powder Pillows and Sodium Thiosulfate Standard Solution for approximately 100 tests.

---

### Specifikationer

|                    |                           |
|--------------------|---------------------------|
| Mätområde:         | 0.2 - 4 mg/L, 1 - 20 mg/L |
| Metodnamn:         | Iodometric Titration      |
| Specialfunktioner: | CH-65 Test Kit            |