



Reagenssats, klor (total), HR digital titrator

Produktnr: 2444800

Farligt

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.



For determination of HR Total Chlorine by Digital Titration with Sodium Thiosulfate. Range: 2000 to 70,000 mg/L Cl₂. Method: Iodometric Set includes Dissolved Oxygen 3 Powder Pillows, Potassium Iodide Powder Pillows, Starch Indicator Solution and Sodium Thiosulfate Digital Titrator Cartridges (2) 2.00 N. Approximately 100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikationer

Antal tester :	100
Mätområde:	2000 - 70000 mg/L Cl ₂
Metodnamn:	Iodometric Titration
Metodnummer:	8209
Plattform:	Digital Titrator
Specialfunktioner:	HR Digital Titrator