



Reagenssats, hårdhet (total), burett, LR

Produktnr:

2447600

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



Reagent set for the determination of Total Hardness by the ManVer[®] 2 Buret Titration Method (Hach method 8226). EPA Accepted method. Range: 0-2,500 mg/L as CaCO₃. Contains: Buffer Solution (Hardness 1), ManVer[®] 2 Hardness Indicator Powder Pillows and 0.0200 N EDTA (TitraVer[®]) Titrant for up to approximately 100 tests.

Specifikationer

Antal tester :	100
EPA compliant:	Ja
Mätområde:	0 - 2500 mg/L CaCO ₃
Metodnamn:	EDTA Titration
Metodnummer:	8226
Parameter:	Hardness (Total)
Specialfunktioner:	Reagent set, Buret, LR