



Reagenssats, hårdhet (total), burett, HR

Produktnr:

2447700

Farligt

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.



Reagent set for the determination of Total Hardness by the ManVer® 2 Buret Titration Method. Range: 100-25,000 mg/L as CaCO₃. Contains: Buffer Solution (Hardness 1), ManVer® 2 Hardness Indicator Powder Pillows and 0.200 N EDTA (TitraVer®) Titrant for up to approximately 100 tests.

Specifikationer

Antal tester :	100
EPA compliant:	Ja
Mätområde:	100 - 25000 mg/L CaCO ₃
Metodnamn:	EDTA Titration
Metodnummer:	8222
Parameter:	Hardness (Total)
Specialfunktioner:	Reagent set, Buret, HR