



Reagenssats, klor (fritt och totalt), DPD, 10 mL

Produktnr:

2459500

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



For determination of Chlorine (Free and Total) by the DPD Method. Hach Methods 8021 and 8167. Range: 0.02 - 2.00 mg/L Cl_2 . Sample Size: 10 mL. Reagent Set contains DPD Free Chlorine Reagent and DPD Total Chlorine Reagent Powder Pillows for 100 tests each.

Boric Acid free!

Specifikationer

Antal tester :	100
Mätområde:	10 mL
Metodnamn:	DPD