



Flask, Erlenmeyer, Wide Mouth, 2000 mL, Subdivisions 600 mL

Produktnr:	2489454
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Graduated glass flask with wide mouth for easier titrations. Reinforced tops with round finish provide maximum strength. Durable graduations indicate approx. volume. Stopper not included. Requires Size 13 rubber/silicone stopper.

Specifikationer

Enhet:	1 each
Kapacitet :	2000 mL
Material:	Glass Wide Mouth