



Sample Container, 120 mL High PF

Produktnr:	2552200
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

These 120 mL (4oz.) High Profile Corning Snap-Seal™ polypropylene containers are designed for a wide variety of applications and provide a reliable leak-tight seal when closed properly. Approximately 3.5 inches tall (91 mm) and 1.75 inches (45 mm) diameter. Plastic, snap seal, translucent. Sold individually.

Specifikationer

Diameter, yttre:	45 mm (approx)
Kapacitet :	120 mL - 1 each
Volym/förpackningsstorlek:	1 /pk