



Sample container, 240 mL

Produktnr:	2552411
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Plastic, leak-proof. Corning Snap-Seal offers a non-sterile, leak-proof, attached cap container for field use in sampling water, soils, and sludge. Locking device helps prevent accidental opening of sample; container closes with unique sealing mechanism so sample won't leak.

Specifikationer

Diameter, yttre:	80 mm (approx)
Kapacitet :	240 mL
Volym/förpackningsstorlek:	140 /pk