



Bottle, wash, diamond alert, 500 mL, Variety 6/pk

Produktnr:	2613049
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Wide-mouthed and color-coded wash bottles provide safe and positive identification of common laboratory solutions. Standard Diamond Alert symbols indicate type of hazard contents may represent. Polyethylene bottles with color-coded polypropylene screw caps and spray-tip assembly. Variety 6-pack includes one bottle each for water, acetone, methanol, isopropanol, toluene, and ethanol.

Specifikationer

Färg:	Transparent
Kapacitet :	500 mL
Material:	LDPE
Stem:	Blue Cap/Stem
Unit:	6 pk
Volym/förpackningsstorlek:	6 pk