



Filter, membrane, pore size 0.2 µm, diameter 47 mm, 200/pk (Gelman)

Produktnr:	2618800
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

For isolation and enumeration of stressed organisms (*Pseudomonas* species). High retention and flow rate. Certified for use in membrane filtration analysis. White polyethersulfone media. STerile and individually packaged. Pore size: 0.2 µm. Diameter: 47mm.

Specifikationer

Diameter:	47 mm
Filter Type:	Filter
Material:	Polyethersulfone
Porstorlek:	0.2 µm
Volym/förpackningsstorlek:	200 /pk