



Media set, replacement, MEL/MF field filtration laboratory

Produktnr: 2619101

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Replacement Media Set for MEL/mColiBlue24 Field Filtration Lab (Product # 2586400).

Media set contains m-ColiBlue24 broth in plastic ampules (4 boxes of 50 ampules), Whirl-Pak Sampling Bags with thiosulfate, sterile, disposable petri dishes with pads, push fit funnel assembly and membrane filters.