



Sample vial, 120 mL w/out thiosulfate, 200/pk

Produktnr:	2634002
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene.

Specifikationer

Diameter:	43 mm
Kapacitet :	Fill line at 100 mL ($\pm$ 2.5 mL)
Längd:	93 mm
Specialfunktioner:	Sample Vial, without Thiosulfate
Volym/förpackningsstorlek:	200/pk