



Thermometer, Enviro-Safe, Non-Mercury, -10 - 260 °C

Produktnr:	2635701
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Enviro-Safe non-Mercury thermometers feature an exclusive biodegradable liquid and dye. Easy to read, aesthetically distinctive due to dark green liquid and permanently fused markings. Each thermometer, individually calibrated against NIST traceable equipment, identified by unique serial number for ISO

9000 or other quality systems.

Specifikationer

Längd:	405 mm
Temperaturområde :	-10 °C - 260 °C
Termometer av vätsketyp:	Non-Mercury
Termometerindelningar:	1 °C
Volym/förpackningsstorlek:	1 each