



Thermometer, Enviro-Safe, Non-Mercury, -20 - 110 °C

Produktnr: 2635702
SEK Pris: Kontakta oss
Inget leveransdatum anges

Enviro-Safe non-Mercury thermometers feature an exclusive biodegradable liquid and dye. Easy to read, aesthetically distinctive due to dark green liquid and permanently fused markings. Each thermometer, individually calibrated against NIST traceable equipment, identified by unique serial number for ISO9000 or other quality systems.

Specifikationer

Temperaturområde :	-20 - 110 °C
Termometer av vätsketyp:	Non-Mercury
Termometerindelningar:	1 °C
Volym/förpackningsstorlek:	1 each