



Cylinder, Graduated, 1000 mL, Certified

Produktnr:	2636253
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in white marking. Bumper guard included for all sizes ≥ 25 mL.

Specifikationer

Certifieringar :	Class A
Höjd :	465 mm
Kapacitet :	1000 mL
Material:	Glass
Mätmetod:	To Deliver
Typ:	Single Metric
Volym/förpackningsstorlek:	1 each