



Cylinder, graduated, mixing, certified, 50 mL

Produktnr:	2636341
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

Specifikationer

Certifieringar :	Class A
Kapacitet :	50 mL
Material:	Glass
Mätmetod:	To Contain
Typ:	Single Metric, Mixing
Volym/förpackningsstorlek:	1 each