



Flask, volumetric, class A, certified, 1000 mL

Produktnr: 2636653

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Individually serialized glass flask supplied with Certificate of Graduation Accuracy and traceable to NIST Standards. Class A tolerances meet requirements of ASTM E438 and ASTM E288. Graduation ring on neck of flask. Supplied with ground glass stopper. (Stopper size 22).

Specifikationer

Certifieringar :	Yes
Enhet:	1 each
Kapacitet :	1000 mL
Material:	Class A Glass