



Iron (Ferrous) Color Disc Test Kit, Model IR-18C, 0-7 mg/L

Produktnr:

2667200

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



The most accurate visual test kit method, Hach color discs provide fast, accurate comparisons. Range: 0-7 mg/L. Contains 100 tests.

Simple, fast, color disc method

Rugged test kit contains all required reagents and apparatus to perform analysis in the field

Specifikationer

Antal tester :	100
Kapslingstyp:	B
Mätområde:	High Single Range
Metodnamn:	Color disc/ 1,10 Phenanthroline
Modell:	IR-18C
Parameter:	Iron, medium range - As Fe ²⁺ (ferrous iron)
Plattform:	Color Disc
Smallest increment:	0.2
Test-kit modell:	IR-18C
Typ:	Single Metric