



Test kit, hypochlorite, model CN-HR

Produktnr:

2687200

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Hypochlorite Kit, Model CN-HR - Drop Count

Utilizing the simple, economical drop count titration method, this kit is used to determine high range hypochlorite solution concentrations (bleach solutions). Range: 5-15%, contains 100 tests.

Simple, economical titration method

Kit contains all required reagents, solutions, and glassware

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	50 - 150 g/L Cl_2 (5-15%)
Metodnamn:	Drop count titration/ Thiosulfate
Parameter:	Hypochlorite, high range - (Bleach) As Cl_2
Plattform:	Drop Count Titration
Smallest increment:	0.5%
Test-kit modell:	CN-HR