



Cleaning Solution for APA6000 Nitrate, 1 L

Produktnr:	2697453
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Specifikationer

Instrument:	APA6000
Lagringsförhållanden:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Mätområde:	N/A
Parameter:	Hardness
Volym/förpackningsstorlek:	1000 mL

Obligatoriska tillbehör

- Skyddshandskar storlek L, blå, nitril, puderfria, 50 styck (Item SM743L)