



## Economy kit m-ColiBlue24®, glass ampules, 200 test

<b>Produktnr:</b>	<b>2779200</b>
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

In 24 hours, Hach's m-ColiBlue24® Broth gives you simultaneous results for total coliforms and E. coli. Use m-ColiBlue24® Broth to accurately monitor and evaluate total coliform and E. coli in drinking water; wastewater; bottled water; beverages; surface, ground, and well water; and ultrapure, chemical processing and pharmaceutical processing waters.

With ready-to-use m-ColiBlue24® Broth, you can differentiate colonies easily - red and blue indicate total coliforms and blue specifies E. coli. Get optimal recovery of stressed and injured organisms. Minimize background growth of non-coliform bacteria.

**Read and confirm results in 24 hours**

**Enumerate total coliforms and E. coli on one petri dish**

**Conduct emergency testing and routine monitoring**

**Get superior sensitivity - 1 CFU/100mL**

**Eliminate special equipment and ultraviolet lamp**