



Hardness Indicator Solution, 0.3 mg/L, 500 mL

Produktnr:

2794649

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



For use with the SP510 Hardness analyser.

Specifikationer

Instrument:	SP510
Koncentration:	0.3 mg/L
Lagringsförhållanden:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Parameter:	Hardness

Obligatoriska tillbehör

- Skyddshandskar storlek L, blå, nitril, puderfria, 50 styck (Item SM743L)